

Erba Mannheim

XL SYSTEM PACKS



XL-640



XL-600



XL-300



XL-1000



XL-200



XL-100

Get the best out of your XL Analyzer

XL SYSTEM PACKS

...are bar coded, ready for use reagents that can be directly loaded on to the XL 200/300/600/640 reagent trays.

Sr. No.	PRODUCT CODE	PRODUCT NAME	METHOD	LINEARITY	PACK SIZE
1	XSYS001	ALBUMIN	BCG	6 g/dl	10 X 44 ml
2	XSYS002	ALP	IFCC, AMP	700 U/L	R1-2 X 44 ml ; R2-2 x 11 ml
3	XSYS003	AMYLASE	CNPG3	1500 U/L	5 X 22 ml
4	XSYS028	BID	DIAZO	25 mg/dl	R1-6X 44ml ; R2-3 X 22ml
5	XSYS023	BIT	DIAZO	25 mg/dl	R1-6X 44ml ; R2-3 X 22ml
6	XSYS007	CALCIUM (A)	ARSENAZO	16 mg/dl	10 x 12 ml
7	XSYS008	CHLORIDE	MERCURIC THIOCYANATE	130 mEq/L	10 X 12 ml
8	XSYS009	CHOLESTEROL	CHOD POD	750 mg/dl	10 X 44ml
9	XSYS022	CK	DGKC & IFCC	1000 U/L	R1-2 X 44 ml ; R2-2 x 11 ml
10	XSYS029	CK MB	IMMUNOINHIBITION	1000 U/L	R1-2 X 44 ml ; R2-2 x 11 ml
11	XSYS024	CREATININE	ALKALINE PICRATE	20 mg/dl	R1-5 X 44 ml ; R2-5 X 11ml
12	XSYS011	GGT	SZASZ	500 U/L	R1-2 X 44 ml ; R2-2 x 11 ml
13	XSYS012	GLUCOSE	GOD POD (with mutarotase)	500mg/dl	10 X44 ml
14	XSYS045	HBA1C	Enzymatic	Upto 16%	Lysing Reagent.: 1x40ml Reagent R1a : 1x17ml Reagent R1b: 1x8ml Reagent R2: 1x11ml Blanking Solution: 1ml Calibrator: 2x1ml
15	XSYS043	HDL DIRECT	Modified PEGME	184.8 mg/dl	R1-4 x 30 ml ; R2-4 x 10ml
16	FBCEM0075	HDL/LDL Calibrator	N.A	N.A.	N.A.
17	XSYS013	LDH-P	DGKC	1000 U/L	R1-2 X 44 ml ; R2-2 x 11 ml
18	XSYS044	LDL DIRECT	IMMUNOINHIBITION	250mg/dl	R1-2 X 30ml ; R2-2 X 10ml
19	XSYS040	MAGNESIUM	Calmagite	4 mg/dL	2x44 ml
20	XSYS027	MICROPROTEIN	PYROGALLOL RED	400mg/dl	10 X 12 ml
21	XSYS015	PHOSPHORUS	UV PHOSPHOMOLYBDATE	15mg/dl	10 X12 ml
22	XSYS016	SGOT	MODIFIED IFCC	240U/L	R1-6X 44ml ; R2-3 X 22ml
23	XSYS017	SGPT	MODIFIED IFCC	240U/L	R1-6X 44ml ; R2-3 X 22ml
24	XSYS018	TOTAL PROTEIN	BIURET	16g/dl	10 X 44 ml
25	XSYS041	TRIGLYCERIDES-SR	GPO	600 mg/dL	10 x 44 ml
26	XSYS020	UREA	UREASE GLDH	300mg/dl	R1-5 X 44 ml ; R2-5 X 11ml
27	XSYS042	URIC ACID - SR	URICASE	25 mg/dl	10 x 44 ml
28	XSYS046	ASO*	IMMUNOTURBIDIMETRIC	800 IU/ml	R1:2x40ml ; R2:2x10ml
29	XSYS047	CRP*	IMMUNOTURBIDIMETRIC	150 mg/L	R1:2x40ml ; R2:2x10ml
30	XSYS048	RF*	IMMUNOTURBIDIMETRIC	160 IU/ml	R1:2x40ml ; R2:2x10ml
31	REKT901	ERBA NORM	NA	NA	4 x 5 ml
32	REKT902	ERBA PATH	NA	NA	4 x 5 ml
33	REKT910	ERBA XL WASH	NA	NA	4 x 100 ml
34	XSYS034	XL MULTICAL	NA	NA	4 X 3 ml